

REV. 3-27-06

65 PLAIN AVENUE NEW ROCHELLE, NY 10801 914-235-9100 • FAX: 914-235-9697 800-473-6682 ortho-rite.com

FOR LAB USE ONLY

THE STABILIZER		
OFFICE NAME	DATE	
DOCTOR		
ADDRESS		
PATIENT NAME		
SHOE SIZE TYPE WORN		
DIAGNOSIS AND		
OBSERVATIONS		
CAST TAKENNON WEIGHT BEARING		EARING
BRACE RIGHT LEFT BILATERAL	1 3/16" PPT PLASTAZOTE 1 1/8" SPENCO	
POST SPECIFICATIONS REARFOOT	☐ 3/16" MULTI-COLOR MEDIUM DENSITY EVA TOP COVER	
☐ EXTRINSIC ☐ INTRINSIC ☐ Post according to lab evaluation ☐ 0 degrees ☐ 3 degrees ☐ 5 degrees	□ ADD FIRST RAY CUT OUT	
FOREFOOT EXTRINSIC	FOR AN ADDITIONAL CHARGE	
☐ Post according to lab evaluation	PPT ARCH FILL	
□ 0 degrees □ 3 degrees □ 5 degrees	☐ TEMPORARY 90 DEGREE ANKLE LOCK	
HEEL CUP DEPTH ☐ 10 mm ☐ 15 mm ☐ 35 mm (standard)	☐ ADDITIONAL MEDIAL/LATERAL INTERFACE F☐ PPT GEL COLD THERAPY INSERTS	ADS
REQUIRED MEASUREMENTS		
A) Width of Forefoot at the Metatarsal Heads		
inches B) Width of Ankle Joint at the Widest Point	STANDARD SPECIFICATIONS	IOOEO
inches	CUSTOM MOLDED POLYPROPYLENE FOOT ORTH FREE MOTION AND F. 101NT ATTACHED TO A UNIT	
C) Circumference above the ankle	 FREE MOTION ANKLE JOINT ATTACHED TO A UNI POLYPROPYLENE INTEGRATED ON PIECE CALF I 	
inches	35 MM DEEP HEEL CUP	
D) Circumference at 3" above the ankle	REMOVABLE MEDICAL AND LATERAL INTERFACE	PAD
inches LAB CAST CORRECTION FEE*	BALANCED BIODYNAMIC FOREFOOT AND	
*To avoid additional charges, we suggest casts be	REARFOOT SYSTEM	
taken with the ankle at 90° Neutral and the Foot in Subtalar Neutral	STANDARD MET LENGTH SHELL WITH COMPRESSIBLE POST TO SULCUS	
Special Instructions	VELCRO® CLOSURE	
opecial instructions_	PPT® ULTRALUX® PADDED	
	TOP COVER TO SULCUS	